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
**Experiment 5.41 2-METHYLHEXAN-2-OL** (*tertiary alcohol from a ketone*)

The apparatus and experimental details are similar to those given in the previous sections. Prepare a Grignard reagent from 24.5 g (1 mol) of magnesium turnings, 137 g (107 ml, 1 mol) of butyl bromide and 450 ml of sodium-dried ether. Add slowly with rapid stirring, and cooling with ice if necessary, a solution of 58 g (73.5 ml, 1 mol) of dry acetone in 75 ml of anhydrous ether. Allow the reaction mixture to stand overnight. Decompose the product by pouring it on to 500 g of crushed ice; dissolve the precipitated magnesium compound by the addition of 10 per cent hydrochloric acid or of 15 per cent sulphuric acid. Transfer to a separatory funnel, remove the ether layer and extract the aqueous solution with three 50 ml portions of ether. Dry the combined etheral solutions over anhydrous potassium carbonate or anhydrous calcium sulphate, filter, distil off the ether, and fractionally distil the residue. Collect the 2-methylhexan-2-ol at 137–141 °C. The yield is 105 g (90%). Assign the  $^{13}\text{C}$  n.m.r. absorptions which occur at  $\delta$  14.1, 23.4, 26.8, 29.3, 44.0 and 70.8.